

Determination of *N*-nitroso sertraline in sertraline preparations with LC-MS

1 Objective and scope

The present method is used to quantify *N*-nitroso sertraline in various sertraline preparations by LC-MS/MS.

2 Principle of the method

Various preparations are extracted with MeOH and the *N*-nitroso sertraline content is determined using an external calibration curve. Two preparations are made per sample, with one preparation being spiked with *N*-nitroso sertraline. *N*-nitroso sertraline exists in 2 isomers, which are, however, difficult to separate. For calibration, the areas of both peaks are integrated together and plotted against the concentration.

3 Validation data

See *in situ* validation.

4 Context

- Method development A-027876

5 Definitions and abbreviations

See **OMCL Glossary**

Term / Abbreviation	Description	Additional information
MRM	Multiple Reaction Monitoring	Selective and sensitive MS measuring mode

6 Additional documents

For specific instructions, see **BPM / Procedures / Tools** or **LIMS / Methods**

None.

7 Special measures / Safety instructions

N-nitroso sertraline is potentially carcinogenic. Appropriate protective measures must be taken.

8 Reference and control substances, test equipment, materials, chemicals and solutions

8.1 Reference material

Designation	Content / purity	LIMS S-No.	Manufacturer/supplier/art. no. (e.g.)
<i>N</i> -nitroso sertraline	0.9983	S-4669-001	SynZeal Research Private Limited (SZ-S002026)

8.2 Control substances

Not applicable.

8.3 Equipment and materials

Designation	LIMS procedure documentation E: documentation of results S: documentation of substances used
LC-MS/MS (e.g. MS (AB Sciex QTRAP 5500) / 0524A	E

8.4 Chemicals

Name	S-No. LIMS	Manufacturer/supplier/art. no. (e.g.)
Acetonitrile	S-2058	AppliChem GmbH / Axon Lab AG / 221881.1611
Methanol	S-1712	AppliChem GmbH / Axon Lab AG / 221091.1612
Formic acid	S-2031	Merck / Merck / 253
Milli-Q Water	S-2206	OMCL

8.5 Solutions

Solutions	Preparation	Shelf life / storage temp.
<i>N</i> -nitroso sertraline stock solution	Transfer about 1 mg of <i>N</i> -nitroso sertraline into a 1.5 mL HPLC vial and dilute to 1 mg/mL with MeOH	7 days at -20°C
<i>N</i> -nitroso sertraline Dil 0	Dilute 100 µL of the stock solution to 10 mL MeOH in a volumetric flask	2 days at 2-8°C
<i>N</i> -nitroso sertraline Dil 1	Dilute 100 µL of the Dil 0 solution to 10 mL MeOH in a volumetric flask	2 days at 2-8°C

To rule out errors in the weighing or dilution of the solution, the stock solution and diluted solutions (Dil 0 and Dil 1) are prepared in duplicate (A/B). As a control, a solution of the calibration curve (approx. 10 ng/mL level) from preparation series B is prepared in duplicate and compared with the calibration curve from preparation series A. The calculated concentration of the sample must not deviate by more than $\pm 5\%$ from the target concentration.

9 Procedure

A mixed sample is prepared from 6 tablets. Weigh out a quantity of the mixed sample that corresponds to 5 mg of active ingredient. However, the maximum quantity of medicinal product to be used is 150 mg. Two preparations are made per sample, one of which is spiked with *N*-nitroso sertraline (approx. 5 ng/mL). The 2 samples are extracted with MeOH for 5 minutes and then centrifuged and the supernatant taken up (with a pipette). If a measured value is outside the calibration range, the sample mass is either increased or decreased. The initial mass is changed so that the expected measured value is within the calibration curve. Alternatively, the measured solution can be brought into the calibration range by dilution.

For quantification, a linear calibration curve with 6 points is produced over the range from approx. 1 ng/mL to approx. 100 ng/mL (approx. 0.2 - 20 ppm). The peak area is plotted against the concentration.

The reference standard (calibration solution 4) is re-analysed at regular intervals to check the instrument performance.

LC system

Column	Luna 3u Phenyl-Hexyl, 50x2 mm, 3 µm		
Mobile phase A	Milli-Q water with 0.1% formic acid		
Mobile phase B	Acetonitrile with 0.1% formic acid		
Autosampler temperature	5°C		
Column temperature	40°C		
Injection volume	8 µL		
Flow rate	0.3 mL/min		
Gradient	RT/ min	%A	%B
	0	40	60
	1.0	40	60
	4.0	40	60
	4.1	0	100
	4.6	0	100
	4.7	40	60
	7.0	40	60
UV Detection (optional)	220 nm, 4 nm bandwidth		

MS settings

Source	ESI / positive
Scan type	MRM
MRM detection window	Unscheduled
Curtain gas (CUR)	35
Collision Gas (CAD)	Medium
Temperature (TEM)	350
Ion Source Gas 1 (GS1)	35
Ion Source Gas 2 (GS2)	10
Ion Spray Voltage [IS]	5500
Entrance Potential (EP)	10
Collision Cell Exit Potential (CXP)	13
Declustering Potential (DP)	60

MRM transitions

ID	Q1	Q3	CE	Dwell time [msec]
NA-Sertraline _Quant	335.1	275.1	11	200
NA-Sertraline _Qual1	335.1	159.0	34	200
NA-Sertraline _Qual2	335.1	129.1	30	200

N-nitroso sertraline elutes at the retention time of 3.1 min

Diverter valve

RT	Diverter Valve
0.0	to waste
1.5	to MS
6.0	to waste

10 Evaluation and measurement uncertainty

10.1 Evaluation

To determine the *N*-nitroso sertraline content, the area of the unspiked sample is converted into the concentration using the calibration line (6-point).

The following formulae are used to calculate the recovery:

$$recovery \% = \frac{C_{sp} - \frac{C_0 \cdot M_{sp}}{M_0}}{C_s} \cdot 100$$

- C_{sp} : calculated concentration of the spiked sample / ng/mL
- C_0 : calculated concentration of the unspiked sample / ng/mL
- M_{sp} : mass of the spiked sample / mg
- M_0 : mass of the unspiked sample / mg
- C_s : concentration of the added spike in the sample / ng/mL

10.2 Measurement uncertainty

Calculated as part of the *in-situ* validation.

11 Documentation and protocols

The requirements for documentation are given in the corresponding operating instructions on this topic. If necessary specific requirements for documentation of process variable parameters are listed below.

No specific documentation is required for process variables.

12 Quality control

For quality control purposes, SST criteria are defined for the test series as well as for each sample.

SST criteria for analytical test series

Parameters	Criterion
Linearity Blank spikes (start of the sequence)	$r \geq 0.995$
Recovery reference standard	70% - 130%

SST criteria for Sample

Parameters	Criterion
Spike recovery	70% - 130%
Qualifier ratios	$\pm 25\%$

13 Document history

The last 5 version changes are shown on the document; previous versions must be deleted from the table.

Version no:	Modification date/approval:	Change from previous version:
01		Creation