
User Instructions for Respirator



With the valve on the left open, place mask over face to check for the right fit.



Fit Upper Strap Stretch the upper strap over your head, above the ear.



Adjust Lower Strap Pull the lower strap and fit it over your head, below your ear.



Press on nose clip with both hands to ensure tight fit around face.



Blow To Check If air leaks from edges, readjust mask and straps to correct fit.



1. Before using, check the respirator for holes, tears or other damage to its surface, strap and/or nose cushioning to ensure that your respirator works well. Replace with a new respirator, if damaged or damp.

2. Select the right respirator size to fit your face. Consult the given size chart on the box for size selection. If your selected size fails to fit your face properly, choose another size until the proper size to fit your face is obtained. Please refer to procedures 3 to 5 below.

3. Fold the nose clip into a _____ shape for a proper fit.

4. The respirator strap can be adjusted for proper fit. Pull the strap outwards at either side of the buckle for tighter fit; push the strap inwards at either side of the buckle for looser fit.

5. To ensure a tight seal on your face, cover the middle surface of the respirator with one hand and exhale deeply. Do not push the respirator against your face. If you feel air blowing at your face or eyes, the respirator's position and strap needs to be re-adjusted. Repeat the above steps until you do not feel air blowing at your face or eyes.

6. For use of the respirator with the micro ventilator, Fit check should be performed with the Micro Ventilator attached.