

## 1-primer mutagenesis

### Primer design

Design one forward primer to anneal at 55°C both upstream and downstream of the point mutation (approximately 45bp total primer length). Use the primer design spreadsheet at <http://macrolab.berkeley.edu> if you need it. Change the fewest number of bases possible to get the desired amino acid.

### Mutagenesis reaction

400ng plasmid (~2uL)  
2uL of 5uM primer  
0.4uL 100mM dNTPs  
1uL PfuTurbo  
5uL Pfu buffer  
ddH<sub>2</sub>O to 50uL total volume

Note: if you order your oligo from a commercial source, it was already run over a desalting column. No further purification is necessary.

95° – 3min  
    95° – 30s  
    53° – 1min                   18 cycles  
    68° – 2min/kb  
68° – 4min/kb

### DpnI digest

4uL mutagenesis reaction  
2uL NEB buffer 4  
0.5uL DpnI  
13.5uL ddH<sub>2</sub>O

Incubate at 37C for 2-2.5 hours

Transform cloning strain of E. coli (e.g. XL-1 Blues) with 6uL of DpnI rxn (30min on ice, 45s heat shock at 42C, 2 min on ice, add 500uL 2YT and plate 200uL on appropriate antibiotic plate).

I get ~400cols/plate, ~45% mutants and 55% wild type plasmids using this protocol.

Things NOT to do:

Don't use Herculase II

Don't undercut the extension time

Don't use any buffer other than the commercially supplied one